


Stone Extraction Baskets

Single Use

Stone Extraction Baskets

- 45° Injection Port & Integrated Sealing Cap
User-friendly and Easy Injection of Contrast Medium
- Detachable Instrument Holder & Outer Sheath
Applicable to Emergency Lithotripsy Spiral and StoneBuster Handle in Emergency Cases
- Ergonomic Instrument Holder
Simple and Safe Basket Manipulation



Stone Extraction Baskets

Single Use

Features & Benefits

- 45° Injection Port & Integrated Sealing Cap
User-friendly and Easy Injection of Contrast Medium
- Detachable Instrument Holder & Outer Sheath
Applicable to Emergency Lithotripsy Spiral and StoneBuster Handle in Emergency Cases
- Ergonomic Instrument Holder
Simple and Safe Basket Manipulation

Stone Extraction Basket Set – Set Consists of Basket, Outer Sheath (Working Length 180 cm), Sealing Cap and Instrument Holder

Catalog No.	Description	Basket Diameter	Basket Length	Outer Sheath Diameter
GEB-11-19-420	4-Wire Baskets	20 mm	42 mm	1.9 mm
GEB-11-19-425	4-Wire Baskets	25 mm	52 mm	1.9 mm
GEB-11-19-430	4-Wire Baskets	30 mm	61 mm	1.9 mm
GEB-11-19-620	6-Wire Baskets	20 mm	42 mm	1.9 mm
GEB-11-19-625	6-Wire Baskets	25 mm	52 mm	1.9 mm
GEB-11-19-630	6-Wire Baskets	30 mm	61 mm	1.9 mm
GEB-11-25-420	4-Wire Baskets	20 mm	42 mm	2.5 mm
GEB-11-25-425	4-Wire Baskets	25 mm	52 mm	2.5 mm
GEB-11-25-430	4-Wire Baskets	30 mm	61 mm	2.5 mm
GEB-11-25-620	6-Wire Baskets	20 mm	42 mm	2.5 mm
GEB-11-25-625	6-Wire Baskets	25 mm	52 mm	2.5 mm
GEB-11-25-630	6-Wire Baskets	30 mm	61 mm	2.5 mm

Distributed by:

Manufactured by:

Medi-Globe GmbH | Medi-Globe-Straße 1-5 | 83101 Achenmühle, Germany
☎ +49 80 32 / 973- 379 | 📠 +49 80 32 / 973- 339 | ✉ sales@medi-globe.de