

# 8-Wire Nitinol Stone Extraction Basket

Single Use

- Kink Resistant Nitinol Baskets  
*Reliable and controlled radial deployment in the bile duct*
- Rotatable Basket in 8-Wire Design  
*Synchronous basket <-> Handle rotation*
- Secure capture and optimal retrieval of small stones and fragments*
- Double Lumen Guidewire Compatible Versions  
*Provide excellent ductal access*
- Precise maneuverability & positioning in challenging bile duct anatomy*



MEDI-GLOBE  
GROUP

**Medi-Globe®**

*The Spirit of Care.*

# 8-Wire Nitinol Stone Extraction Basket

Single Use

## Features & Benefits

- Kink Resistant Nitinol Baskets  
*Reliable and controlled radial deployment in the bile duct*
- Rotatable Basket in 8-Wire Design  
*Synchronous basket <> Handle rotation  
Secure capture and optimal retrieval of small stones and fragments*
- Double Lumen Guidewire Compatible Versions  
*Provide excellent ductal access  
Precise maneuverability & positioning in challenging bile duct anatomy*

## Double Lumen Guide Wire Compatible Stone Extraction Basket

### 8-Wire Nitinol Stone Extraction Basket – Working Length 180 cm

Article-No.	Description	Basket Diameter	Max. Accepted Guidewire	Outer Diameter	Min. Working Channel
GEB-14-27-820	8-Wire Basket	20 mm	.035"	2.7 mm	3.2 mm
GEB-14-27-830	8-Wire Basket	30 mm	.035"	2.7 mm	3.2 mm

## Single Lumen Stone Extraction Basket

### 8-Wire Nitinol Stone Extraction Basket – Working Length 180 cm

Article-No.	Description	Basket Diameter	Max. Accepted Guidewire	Outer Diameter	Min. Working Channel
GEB-13-19-820	8-Wire Basket	20 mm	-	1.9 mm	2.8 mm
GEB-13-19-830	8-Wire Basket	30 mm	-	1.9 mm	2.8 mm

Caution: This device is not detachable and must not be used with the StoneBuster® Lithotripter or any other mechanical lithotripsy system.

Distributed by:

L-6EB-SU-02-ENG Ver. 2.0 (2017-04)