

ECO-Injectra

Injection Needles – Single Use

ECO-Injectra

- **Metal Tip Protection**
Minimizes the risk of sheath puncture and secures a safe needle guidance in angulated endoscope positions
- **Transparent Needle Tube**
Excellent flow control of medicine and fluids
- **Ergonomic Handle & Needle Lock Mechanism**
Stable & controlled puncture in targeted tissue and vessels

ECO-Injectra

Injection Needles – Single Use

Bronchoscopy

ECO-Injectra – Injection needle

Catalog No.	Needle Diameter	Needle Length	Sheath Outer Diameter	Working Length
GSN-21-18-120	25 Gauge	4 mm	1.8 mm	120 cm
GSN-22-18-120	25 Gauge	6 mm	1.8 mm	120 cm

Gastroscopy

ECO-Injectra – Injection needle

Catalog No.	Needle Diameter	Needle Length	Sheath Outer Diameter	Working Length
GSN-21-23-180	25 Gauge	4 mm	2.3 mm	180 cm
GSN-22-23-180	25 Gauge	6 mm	2.3 mm	180 cm
GSN-23-23-180	22 Gauge	4 mm	2.3 mm	180 cm
GSN-24-23-180	22 Gauge	6 mm	2.3 mm	180 cm
GSN-25-23-180	19 Gauge	4 mm	2.3 mm	180 cm
GSN-26-23-180	19 Gauge	6 mm	2.3 mm	180 cm

Colonoscopy

ECO-Injectra – Injection needle

Catalog No.	Needle Diameter	Needle Length	Sheath Outer Diameter	Working Length
GSN-21-23-230	25 Gauge	4 mm	2.3 mm	230 cm
GSN-22-23-230	25 Gauge	6 mm	2.3 mm	230 cm
GSN-23-23-230	22 Gauge	4 mm	2.3 mm	230 cm
GSN-24-23-230	22 Gauge	6 mm	2.3 mm	230 cm
GSN-25-23-230	19 Gauge	4 mm	2.3 mm	230 cm
GSN-26-23-230	19 Gauge	6 mm	2.3 mm	230 cm