

Polyp & Foreign Body Retrievers

Single Use

Polyp & Foreign Body Retrievers

- Strong & Wide Opening Graspers
Large Variety of Applications
- Ergonomically Shaped Handle
User-Friendly and Safe




 **Medi-Globe**[®]
The Spirit of Care.[™]

Polyp & Foreign Body Retrievers

Single Use

Polyp & Foreign Body Retrievers



Foreign Body Retriever – Grasper with Hooks

| Catalog No. | Grasper Shape | Catheter Outer Diameter | Working Length |
|---------------|-----------------|-------------------------|----------------|
| GFB-01-25-180 | 4-Prong Grasper | 2.5 mm | 180 cm |
| GFB-01-25-230 | 4-Prong Grasper | 2.5 mm | 230 cm |

Polyp Retriever – Grasper with Eyelets

| Catalog No. | Grasper Shape | Catheter Outer Diameter | Working Length |
|---------------|-----------------|-------------------------|----------------|
| GFB-06-25-180 | 3-Prong Grasper | 2.5 mm | 180 cm |
| GFB-06-25-230 | 3-Prong Grasper | 2.5 mm | 230 cm |
| GFB-02-25-180 | 4-Prong Grasper | 2.5 mm | 180 cm |
| GFB-02-25-230 | 4-Prong Grasper | 2.5 mm | 230 cm |

Foreign Body Retriever – Grasper with Forked Jaws

| Catalog No. | Grasper Shape | Catheter Outer Diameter | Working Length |
|---------------|-----------------|-------------------------|----------------|
| GFB-03-25-180 | Forked Jaws 2:2 | 2.5 mm | 180 cm |
| GFB-03-25-230 | Forked Jaws 2:2 | 2.5 mm | 230 cm |
| GFB-04-25-180 | Forked Jaws 1:2 | 2.5 mm | 180 cm |
| GFB-04-25-230 | Forked Jaws 1:2 | 2.5 mm | 230 cm |

Foreign Body Basket

| Catalog No. | Grasper Shape | Basket Diameter | Catheter Outer Diameter | Catheter Outer Diameter |
|---------------|---------------|-----------------|-------------------------|-------------------------|
| GFB-07-25-425 | 4-Wire Basket | 25 mm | 2.5 mm | 230 cm |
| GFB-07-25-625 | 6-Wire Basket | 25 mm | 2.5 mm | 230 cm |

Foreign Body Forceps with coating

| Catalog No. | Grasper Shape | Catheter Outer Diameter | Catheter Outer Diameter |
|---------------|---------------|-------------------------|-------------------------|
| GFB-09-23-120 | Rat Tooth | 2.3 mm | 120 cm |
| GFB-09-23-180 | Rat Tooth | 2.3 mm | 180 cm |
| GFB-09-23-230 | Rat Tooth | 2.3 mm | 230 cm |

Foreign Body Forceps with coating

| Catalog No. | Grasper Shape | Catheter Outer Diameter | Catheter Outer Diameter |
|---------------|---------------|-------------------------|-------------------------|
| GFB-10-23-180 | Alligator | 2.3 mm | 180 cm |
| GFB-10-23-230 | Alligator | 2.3 mm | 230 cm |